



ALBAcyte® Expanded Rh Negative Antibody Screen (Z464U)



4 x 5 mL

The ALBAcyte® Expanded Rh Negative Antibody Screen contains three RhD negative cells (r'r, r''r, rr) that can detect all required alloantibodies beyond RhD and a fourth cell (presumptive R₀R₀) for detecting anti-D. The screening cell set meets the FDA CFR requirements for antibody detection.

Who should use this kit?

Any hospital delivering babies

- Detect alloantibodies in RhD negative pregnant women who have received Rh Immune Globulin prophylactically.
- One-step process for performing antibody detection in RhD negative women.
- Reduction in the opportunity for human error as all cells are present in one set.

Donor centers and hospitals with sickle cell clinics

- Support phenotyping for Rh and Kell.
- Heterozygous control cells for C, c, E, e and K antigen testing available within the three RhD negative cells.

Laboratories that perform complex antibody identification workups

- Additional source of red cells for user-selected panels.
- Homozygous antigen expression in clinically significant blood groups, i.e. Fy(a+b-), Fy(a-b+), Jk(a+b-), Jk(a-b+), S+s-, S-s+.



TO FIND OUT MORE, PLEASE SCAN THIS CODE

Visit us at alivedx.com or email us at Customer.ServiceUS@alivedx.com



©Alba Bioscience Ltd, 2023. ALBAcyte® is a trademark or registered trademark of Quotient Holdings Financial Company Limited or its subsidiaries in various jurisdictions. Menus and capabilities are subject to change. Not all methods may be available in all territories. Subject to regulatory clearance.