ALBAcyte® REAGENT RED BLOOD CELLS FOR USE IN IDENTIFICATION OF UNEXPECTED ANTIBODIES Expanded Rh Negative Antibody Screen - REF Z464U



	Patient/Samp] _							
Name:	me: DOB:			ABO Rh:		Lot No:	V272354		
Diagnosis/History:	Hospital:	DAT:	Anti-IgG:	Anti-C3		Expiry Date:	2024.06.17	U.S. License 1807	

	Rh-hr				Kell						D	uffy	ŀ	Kidd		Lewis		MNS			Р	Lutheran		Additional Antigens					TEST RESULTS										
Cell #	Rh-hr	Donor	D	С	Ε	С	е	f	٧	C_{m}	K	k	Kp ^a	Kp^b	*Jsa	*Js ^b	Fy	Fy	Jk	a Jk) Le	e ^a Le ^b	М	N	S	s	P1	Lu ^a	Lu ^b	Xg ^a	Wrª	Sp	ecial	Types					Cell #
1	r'r	4154107	0	+	0	+	+	+	NT	0	0	+	0	+	0	+	+	0	0	+	7	. 0	+	+	+	+	+	0	+	+	0								1
2	r''r	2TJVB1E	0	0	+	+	+	+	NT	0	0	+	0	+	NT	+	0	+	+	0	C	+	+	+	0	+	+	0	+	*+	0								2
3	rr	2AA5GBG	0	0	0	+	+	+	NT	0	+	+	0	+	NT	+	0	+	+	0	C	+	+	0	+	0	+	0	+	*+	0								3
4	R_0R_0	1124-947-1	+	0	0	+	+	+	NT	0	0	+	0	+	0	+	0	0	+	0	C	+	+	0	0	+	+	0	+	*+	0								4
	Pati	ent Cells																																					

Notes:

- 1. All cells are DAT negative.
- 2. NT = Not tested.
- 3. The f antigen status has been determined presumptively.
- 4. Cell #4 is designated as R_0R_0 based on the statistical probability for the source donor population but a genotype of R_0 r cannot be excluded.
- * Indicates those antigens whose presence or absence may have been determined using only a single example of a specific antibody.

Conclusions/Further Testing Required:		
Signature/Initials:	Date:	

Product Technical Support Tel: 1-888-228-1990 Web: www.quotientbd.com