



A F	0	DLE		ITIC	
A	3U/	KIIL	J AN	เมอ	<b>ERA</b>

PRODUCT NAME	CELL LINE	UNIT SIZE	CODE	ADDITIONAL INFORMATION
ALBAclone® Anti-A	LA2	10 x 10 mL	Z001	
ALBAclone® Anti-B	LB2	10 x 10 mL	Z011	Does not detect Acquired B
ALBAclone® Anti-A,B	LA2, LB2, ES15	10 x 10 mL	Z021	Detects Ax cells
ALBAclone® Anti-D alpha (IgM)	LDM1	10 x 10 mL	Z031	Does not detect DVI
ALBAclone® Anti-D optimum (IgM/IgM)	LDM1, ESD1M	10 x 10 mL	Z039	Detects DVI by direct test
ALBAclone® Anti-D blend (IgM/IgG)	LDM3, ESD1	10 x 10 mL	Z041	Detects DVI by IAT
ALBAclone® Anti-D reagent control	N/A	10 x 10 mL	Z271	For use with Z031, Z039, Z041
Al BAlect® Anti-A lectin	N/A	3 x 5 ml	7241	Dolichos hiflorus

#### **REAGENT RED CELLS**

> REVERSE GROUPIN	>	RE\	/ERSE	<b>GRO</b>	UPIN
-------------------	---	-----	-------	------------	------

DPODLICT NAME

PRODUCT NAME	UNIT SIZE	CODE	ADDITIONAL INFORMATION
ALBAcyte® A <sub>1</sub> Cells	1 x 10 mL	Z401	Reverse Grouping & Controls
ALBAcyte® A <sub>2</sub> Cells	1 x 10 mL	Z406	Reverse Grouping & Controls
ALBAcyte® B Cells	1 x 10 mL	Z411	Reverse Grouping & Controls
ALBAcyte® O R <sub>1</sub> r Cells	1 x 10 mL	Z416	Controls
ALBACYTE® Orr Cells	1 v 10 ml	7/121	Controls

## > ANTIBODY DETECTION & IDENTIFICATION

PRODUCT NAME	UNIT SIZE	CODE	ADDITIONAL INFORMATION
ALBAcyte® Antibody Screen (native)	3 x 10 mL	Z451	Patient and Donor Screening
ALBAcyte® RhD Negative CAT Antibody Screen (native)	3 x 5 mL	Z466	Patient Screening Post Prophylactic Anti-D
ALBAcyte® Antibody Identification (native)	10 x 5 mL	Z471	Antibody Identification
ALBAcyte® Antibody Identification (papainised)	10 x 5 mL	Z472	Antibody Identification

ADDITIONAL INFORMATION

LIMIT CIZE

> QUALITY CONTROL RED CELLS			
PRODUCT NAME	UNIT SIZE	CODE	ADDITIONAL INFORMATION
ALBAcyte® IgG Sensitised Cells	1 x 10 mL	Z441	Confirmation & Controls

## QUALITY CONTROL BLOOD GROUP SEROLOGY

•			
PRODUCT NAME	UNIT SIZE	CODE	ADDITIONAL INFORMATION
ALBAsure® Sensitivity Control Kit	4 x 5 mL	Z250	Sensitivity control kit for column agglutination technology – Anti-K, Anti-Ē, Anti-E, Anti-S
ALBAcheck® High Titre Controls Kit	4 x 6 mL	Z257	Positive and negative controls for high titre anti-A/B testing
ALBAcheck® Anti-D <0.1 IU/mL	3 x 5 mL	Z262	Sensitivity control for column agglutination technology and tube testing
ALBAcheck® Anti-D 0.05 IU/mL	3 x 5 mL	Z263	Sensitivity control for column agglutination technology
ALBAcheck® AB Serum (Human)	10 x 5 mL	Z281	Inert Serum – Negative control
ALBAcheck® Simulated Whole Blood Controls	4 x 6 mL	Z489	Controls for manual and automated immunohaematology test platforms

# Building on over 30 years in transfusion diagnostics, Alba supplies thousands of customers around the globe.

## ANTISERA FOR EXTENDED PHENOTYPING

PRODUCT NAME	CELL LINE	UNIT SIZE	CODE	ADDITIONAL INFORMATION*
ALBAclone® Anti-E	DEM1	3 x 5 mL	Z073	Direct Agglutination – 5/15 minutes
ALBAclone® Anti-c̄	H48	3 x 5 mL	Z083	Direct Agglutination – 5/15 minutes
ALBAclone® Anti-k (cellano)	Lk1	1 x 2 mL	Z137	Direct Agglutination – 5 minutes
ALBAclone® Anti- Fy <sup>a</sup>	DG-FYA-02	1 x 5 mL	Z152	IAT – 15 minutes
ALBAclone® Anti- Fy <sup>b</sup>	SpA264LBg1	1 x 5 mL	Z154	Direct Agglutination – Immediate Spin (Detects Fy <sup>x</sup> )
ALBAclone® Anti-Jk <sup>a</sup>	P3HT7	1 x 2 mL	Z162	Direct Agglutination – Immediate Spin
ALBAclone® Anti-Jkb	P3.143	1 x 2 mL	Z166	Direct Agglutination – Immediate Spin
ALBAclone® Anti-M	LM1	1 x 5 mL	Z171	Direct Agglutination – 5 minutes
ALBAclone® Anti-N	LN3	1 x 5 mL	Z176	Direct Agglutination – 5 minutes
ALBAclone® Anti-S	P3S13JS123	1 x 2 mL	Z182	IAT – 10 minutes
ALBAclone® Anti-P1	650	1 x 5 mL	Z202	Direct Agglutination – Immediate Spin
ALBAclone® Anti-Le°	LEA2	1 x 5 mL	Z212	Direct Agglutination – 15 minutes
ALBAclone® Anti-Le <sup>b</sup>	LEB2	1 x 5 mL	Z217	Direct Agglutination – 15 minutes
ALBAclone® Anti-Lu <sup>b</sup>	LU2	1 x 2 mL	Z223	Direct Agglutination – 15 minutes
ALBAsera® Anti-š	N/A	1 x 2 mL	Z186	IAT – 15 minutes
ALBAsera® Anti-Wra	N/A	1 x 2 mL	Z231	IAT – 15 minutes

# \*For reference only – refer to IFU for appropriate technique.

# ANTI-HUMAN GLOBULIN

PRODUCT NAME	UNIT SIZE	CODE	ADDITIONAL INFORMATION
Polyspecific AHG	10 x 10 mL	Z350	Monoclonal Anti-Human C3d / polyclonal rabbit Anti-Human IgG blend (dyed green)
Monospecific Anti-Human IgG	3 x 10 mL	Z356	Polyclonal rabbit Anti-Human IgG (clear)
Monospecific Anti-Human C3d	1 x 3 mL	Z360	Monoclonal Anti-Human C3d (clear)

#### SPECIALTY PRODUCTS

PRODUCT NAME	UNIT SIZE	CODE	ADDITIONAL INFORMATION
ALBAcheck® Competency Testing Kit	10 x 2 mL	Z276	Self-assessment of blood banking techniques and test systems
ALBA Elution Kit	10 elutions	Z313	For use in acid elution of antibodies from intact red blood cells







# **SUPPORTING CUSTOMERS**



Expert technical support from experienced blood bankers



Samples for reagent evaluations available\*

\*Not all products may be available in all territories.



Reliable customer service



TO FIND OUT MORE, PLEASE SCAN THIS CODE

Visit us at alivedx.com or email us at Customer.ServiceEU@alivedx.com

in X

# **Alive**Dx

At AliveDx, we empower diagnostic insights, transform patient care, and innovate for life. With over 30-years in in-vitro diagnostics, we put the health of patients first by creating innovative solutions designed for faster disease detection to shorten the time to diagnosis. Alba—our blood bank reagents portfolio—and MosaiQ®—one of the first multiplexing and multimodality automated testing platforms—are designed to make labs more efficient and clinical decision-making more effective. Our journey started as Alba Bioscience, and later, Quotient. Today, as AliveDx, we innovate for life.

AliveDx's operations are based in Switzerland, Scotland, US and the UAE.

©Alba Bioscience Ltd, 2023. ALBAclone,® ALBAcyte,® ALBAlect,™ ALBAsure,® ALBAhance,™ ALBAzyme,® and ALBAcheck® are trademarks or registered trademarks of Quotient Holdings Financial Company Limited or its subsidiaries in various jurisdictions. Menus and capabilities are subject to change. Not all methods may be available in all territories. Subject to regulatory clearance.