## **ALBAsure**® QC Kit Worksheet



Expected Test Results	Blood Grouping Reagents						Reagent Red Blood Cells										QC Kit Lot Number:		
					Control	۸1		D			Screening Cell 1   Screening Cell 2   Screening Cell 3								
	Anti-A	Anti-B	Anti-AB	Anti-D	Control	A1	A2	В	0	37 °C	IAT	37 ℃	IAT	37 ℃	IAT				
QC Cell 1 – AB rr	1-4+	1-4+	1-4+	0	0											Expiration Date			
QC Cell 2 – O R1r	0	0	0	1-4+	0											1			
QC Antibody						1-4+	1-4+	1-4+	0	0-2+	1-4+	0-2+	1-4+	0-2+	1-4+				
Daily Test	Blood Grouping Reagents								Reagent Red Blood Cells					Additional Reagents					
	Anti-A	Anti-B	Anti-AB	Anti-D	Control	A1	A2	В	0	Screenii 37 °C	ng Cell 1 IAT	Screenii 37 °C	ng Cell 2 IAT	Screeni 37 °C	ng Cell 3 IAT	AHG	Enhancement Media	IgG Sensitize	
Lot:																Lot:	Lot:	Lot:	
Expiration Date:																			
QC Cell 1 – AB rr																Exp:	Exp:	Exp:	
QC Cell 2 – O R1r																			
QC Antibody																			
Satisfactory ✓ Unsatisfactory×																Tech:	Date:		
Lot:																Lot:	Lot:	Lot:	
Expiration Date:																			
QC Cell 1 – AB rr																Exp:	Exp:	Exp:	
QC Cell 2 – O R1r																			
QC Antibody																			
Satisfactory ✓ Unsatisfactory×																Tech:	Date:		
Lot:																Lot:	Lot:	Lot:	
Expiration Date:																			
QC Cell 1 – AB rr																Ехр:	Exp:	Exp:	
QC Cell 2 – O R1r																			
QC Antibody																			
Satisfactory ✓ Unsatisfactory ×																Tech:	Date:		
Lot:																Lot #:	Lot #:	Lot #:	
Expiration Date:																			
QC Cell 1 – AB rr																Ехр:	Exp:	Exp:	
QC Cell 2 – O R1r																			
QC Antibody																			
Satisfactory ✓ Unsatisfactory×																Tech:	Date:		