## ALBAcyte® REAGENT RED BLOOD CELLS FOR DETECTION OF UNEXPECTED ANTIBODIES Antibody Screen (2-Cell) - REF Z454U



Patient/Sample Details										
Name: DOB: Sample Identification: ABO Rh:										
Diagnosis/History:	Hospital:	DAT:	Anti-IgG:		Anti-C3					

Lot No: V283361

Expiry Date: 2025.06.02 U.S. License 1807

			Rh-hr				Kell						Duffy			Kidd		Lewis		MNS			Р	P Lutheran		Additional Antigens				gens	TEST RESULTS								
Cell #	Rh-hr	Donor	D	С	Е	С	е	f	٧	$C_{w}$	K	k	$Kp^{a}$	$Kp^{b}$	*Js <sup>a</sup>	*Js <sup>b</sup>	Fy	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Le	Le <sup>b</sup>	М	N	S	s	P1	Lu <sup>a</sup>	Lu <sup>b</sup>	*Xg <sup>a</sup>	Wr	Spe	ecial	Types					Cell #
1	$R_1R_1$	4487716	+	+	0	0	+	0	NT	0	+	+	0	+	0	+	+	0	+	0	+	0	0	+	0	+	+	0	+	+	0								1
2	$R_2R_2$	2WKZTR*	+	0	+	+	0	0	NT	0	0	+	0	+	NT	+	0	+	0	+	0	+	+	0	+	0	+	0	+	+	0								2
	Patier	nt Cells																																					

## Notes:

- 1. All cells are DAT negative.
- 2. NT = Not tested.
- 3. The f antigen status has been determined presumptively.
- \* Indicates those antigens whose presence or absence may have been determined using only a single example of a specific antibody.

Product Technical Support Tel: 1-888-228-1990 Web: www.alivedx.com

Conclusions/Further Testing Required	d:	
Signature/Initials:	Date:	
Signature/initials.	Date.	

Z454UPI/04