

Why LumiQ®?

Designed to address the needs of laboratories.



Optimized Workflow

Patient sample in, result out, in one single step
True walk-away time



Process Standardization

Minimizes impact of subjectivity and avoids human error
Provides positive and negative pre-classification



Actionable Insights

Image capture and data storage helps patient history management
Provides estimate endpoint titer

//

IFA platform designed to provide fast, actionable insights with an intuitive user experience regardless of level of expertise.

//

Turnkey solution for autoimmune diseases

The combination of MosaiQ® and LumiQ® systems provides a unique, automated, seamlessly integrated solution to cover your needs for screening and identification of autoantibodies for a wide range of diseases.



Find out more about LumiQ®

©AliveDx Suisse SA. 2025. The AliveDx logo, AliveDx, MosaiQ, MosaiQ AiPlex, MosaiQ AllergyPlex COMBO and LumiQ are trademarks or registered trademarks of AliveDx group companies in various jurisdictions. Menus and capabilities are subject to change. Not all methods may be available in all territories. Subject to regulatory clearance. This product is not currently available in the USA, cleared or approved by the U.S. Food and Drug Administration (FDA) for clinical use.

AliveDx

AliveDx



LUMIQ®

All-in-one immunofluorescence solution

Not available for diagnostic purposes in the US. Subject to FDA clearance





Designed to provide workflow efficiency from patient sample to result.

A new level of performance from a benchtop immunofluorescence solution. LumiQ® platform will be able to perform automated IFA testing from patient sample to result.

Indirect immunofluorescence remains the standard method for testing of various antibodies in autoimmune diseases. The detection and measurement of autoantibodies have significantly evolved with automation, offering laboratories more enhanced and standardized options.




Actionable Insights:


-  Provides positive/negative pre-classification
-  Provides estimated endpoint titer
-  Suggests pattern recognition based on sophisticated algorithm
-  Panoramic view of the well

LumiQ® All-in-one solution: Patient sample in, result out.

ANA
HEp-2 cells




dsDNA
Crithidia luciliae




Connective Tissue Diseases

ANCA
Neutrophils Ethanol
Neutrophils Formalin




rLKS
Triple Substrate
Rat Liver/Kidney/Stomach



Autoimmune Vasculitis



Autoimmune Liver Diseases Gastritis



Scan for scientific evidence

One single run, up to:

216
samples

12
different assays

27
slides from different tests

288
dilution capacity

10
patterns recognition