ALBAcyte® REAGENT RED BLOOD CELLS FOR USE IN IDENTIFICATION OF UNEXPECTED ANTIBODIES

Expanded Rh Negative Antibody Screen - REF Z464U



Patient/Sample Details										
Name:	DOB:	Sample Identification	:	ABC	Rh:					
Diagnosis/History:	Hospital:	DAT:	Anti-IgG:		Anti-C3					

Lot No: V288580 Expiry Date: 2025.12.01 U.S. License 1807

			Rh-hr					Kell							Ouffy		Kidd		Lewis		MNS				Р	Lutheran		Additional Antigens					TEST RESULTS								
Cell #	Rh-hr	Donor	D	С	Е	С	е	f	*1	V C)w	K	k	Kp^a	Kp^b	*Js ^a	*Js ^b	, F	^a Fy	Jk	a Jk	b L	e ^a Le	b	М	N	S	s	P1	Lu ^a	Lu ^b	Xg ^a	Wr	^a Sp	ecial	Types					Cell #
1	r'r	426322	0	+	0	+	+	+	N	IT (0	0	+	0	+	0	+	+	0	C	+		- 0		+	0	0	+	+	0	+	+	0								1
2	r"r	1098873	0	0	+	+	+	+	() (0	0	+	0	+	0	+	C	+	+	. 0	() +		+	0	+	0	+	0	+	+	0								2
3	rr	5312271215016	0	0	0	+	+	+	N	IT (0	+	+	0	+	NT	+	C	+	+	+	() +		0	+	0	+	0	0	+	*+	0								3
4	R_0R_0	4426065	+	0	0	+	+	+	N	IT (0	0	+	0	+	0	+	C	0	+	+	() +		+	+	0	+	+	0	+	*+	0								4
	Patie	ent Cells																																							

Notes:

- 1. All cells are DAT negative.
- 2. NT = Not tested.
- 3. The f antigen status has been determined presumptively.
- 4. Cell #4 is designated as R_0R_0 based on the statistical probability for the source donor population but a genotype of R_0r cannot be excluded.
- * Indicates those antigens whose presence or absence may have been determined using only a single example of a specific antibody.

Conclusions/Further Testing Required:		
Signature/Initials:	Date:	

Product Technical Support Tel: 1-888-228-1990 Web: www.alivedx.com

Z464UPI/04